PROBLEMS AND SOLUTIONS (LAB MANUAL, PAGE 26 & 27)

15. Calculate the amount of oxalic acid (C₂H₂O₄.2H₂O, MW 126) for the preparation of 250 ml 0.1 N oxalic acid solutions.

Solution: We know,
$$C = \frac{W \times 1000}{M \times V}$$

Given that, C = 0.1N, V = 250ml, M = 63 (gram-equiv-wt of oxalic acid), W = ?

Therefore, rearranging the above formulae,

$$W = (C \times M \times V)/1000 = (0.1 \times 63 \times 250)/1000 = 1.575 \text{ gm (Ans)}$$

16. Calculate the normality of HCl solution if 5gm of HCl (MW 36.5) is taken in 250 ml solution.

Solution: We know,
$$C = \frac{W \times 1000}{M \times V}$$

Given that, W = 5 gm, V = 250ml, M = 36.5 (gram-equiv-wt of oxalic acid), C = ?

Therefore, putting the values in the above formulae, $C = (5\times1000)/(36.5\times250) = 0.548 \text{ N (Ans)}$

17. If 15ml 0.25 N HCl solutions react with 20ml NaOH solution, then calculate the normality of NaOH solution.

Solution: We know,
$$S_{base} \times V_{base} = S_{acid} \times V_{acid}$$

Given that,
$$V_{acid} = 15 \text{ ml}$$
, $S_{acid} = 0.25 \text{ N}$, $V_{base} = 20 \text{ ml}$, $S_{base} = ?$

Therefore, rearranging the above formulae,

$$S_{base} = (S_{acid} \times V_{acid}) / V_{base} = (0.25 \times 15) / 20 = 0.1875 \text{ N (Ans)}$$

18. Suppose, 9.5 ml of oxalic acid solution and 11.5 ml of dil. HCl solution are required to titrate 10 ml NaOH solution separately. Calculate the normality of NaOH and dil. HCl solutions. [Given that oxalic acid solution is 0.13(N)].

Solution: We know,
$$S_{base} \times V_{base} = S_{acid} \times V_{acid}$$

(A) Given that, $V_{\text{oxalic acid}} = 9.5 \text{ ml}$, $S_{\text{oxalic acid}} = 0.13 \text{ N}$, $V_{\text{NaOH}} = 10 \text{ ml}$, $S_{\text{NaOH}} = ?$

Therefore, rearranging the above formulae,

$$S_{NaOH} = (S_{oxalic acid} \times V_{oxalic acid}) / V_{NaOH} = (0.13 \times 9.5) / 10 = 0.1235 \text{ N (Ans)}$$

(B) Given that,
$$V_{dil.HCl} = 11.5 \text{ ml}$$
, $S_{NaOH} = 0.1235 \text{ N}$, $V_{NaOH} = 10 \text{ ml}$, $S_{dil.HCl} = ?$

Therefore, rearranging the above formulae,

$$S_{dil.HCl} = (S_{NaOH} \times V_{NaOH}) / V_{dil.HCl} = (0.1235 \times 10) / 11.5 = 0.107 N (Ans)$$

19. Calculate the amount of Na₂CO₃ (MW 106) to prepare 0.3 N 200ml Na₂CO₃ solution.

Solution: We know,
$$C = \frac{W \times 1000}{M \times V}$$

Given that, C = 0.3 N, V = 200 ml, M = 53 (gram-equiv-wt of Na₂CO₃), <math>W = ?

Therefore, rearranging the above formulae,

$$W = (C \times M \times V)/1000 = (0.3 \times 53 \times 200)/1000 = 3.18 \text{ gm (Ans)}$$

22. Calculate the amount of Cu^{+2} in a 150 ml blue vitriol solution if 10 ml of it is titrated with 5 ml 0.039N sodium thiosulphate solution. (1 ml 1N $Na_2S_2O_3 = 0.06354$ gm of Cu^{+2}).

<u>Solution</u>: We know, 1 ml 1N Na₂S₂O₃ \equiv 0.06354 gm of Cu⁺²

Given that, V_{thio} = 5 ml, S_{thio} = 0.039 N, volume of blue vitriol solution = 150 ml

Therefore, amount of
$$Cu^{+2}$$
 in 10 ml of blue vitriol solution = $0.06354 \times V_{Thio} \times S_{Thio}$ gm Amount of Cu^{+2} in 150 ml of blue vitriol solution = $0.06354 \times V_{Thio} \times S_{Thio} \times 15$ gm = $0.06354 \times 5 \times 0.039 \times 15$ gm = 0.18585 gm (Ans)

23. Calculate the known value of copper when 3g blue vitriol is dissolved in 100 ml of solution. (Atomic weight of Cu = 63.54, Mol. weight of $CuSO_4.5H_2O = 249.68$).

Solution: We know,
$$\frac{At.wt.of\ Cu\ /\ Fe \times Amount\ of\ salt\ taken\ (in\ gm)}{Mol.wt.of\ blue-vitriol\ /\ Mohr's\ salt}$$

Given that, Atomic weight of Cu = 63.54, Mol. weight of $CuSO_4.5H_2O = 249.68$; Amount of blue vitriol taken = 3 gm;

Therefore, Known value of copper in 100 ml of blue vitriol solution = $(63.54 \times 3)/249.68$ = 0.76345 gm (Ans)

25. Calculate the amount of Fe⁺² in a 300 ml Mohr's salt solution if 10 ml of this solution is titrated with 4 ml 0.075N $K_2Cr_2O_7$ solution. (1 ml 1N $K_2Cr_2O_7 \equiv 0.05584$ gm of Fe⁺²). Ans: 0.50256 gm

Solution: We know, 1 ml 1N $K_2Cr_2O_7 \equiv 0.05584$ gm of Fe⁺²

Given that, $V_{K2Cr2O7} = 4$ ml, $S_{K2Cr2O7} = 0.075$ N, volume of Mohr's salt solution = 300 ml

Therefore, amount of Fe $^{+2}$ in 10 ml of Mohr's salt solution = $0.05584 \times V_{K2Cr2O7} \times S_{K2Cr2O7}$ gm Amount of Fe $^{+2}$ in 300 ml of Mohr's salt solution = $0.05584 \times V_{K2Cr2O7} \times S_{K2Cr2O7} \times 30$ gm = $0.05584 \times 4 \times 0.075 \times 30$ gm = 0.50256 gm (Ans)

--

27. Calculate the amount of Fe^{+2} ions in 500 ml Mohr's salt solution if 10ml of this solution is titrated with 7.5ml 0.06N KMnO₄ solution. (1ml 1N KMnO₄= 0.05584 gm of F^{+2}). Ans: 1.256 gm

Solution: We know, 1 ml 1N KMnO₄ \equiv 0.05584 gm of Fe⁺²

Given that, $V_{KMnO4} = 7.5$ ml, $S_{KMnO4} = 0.06$ N, volume of Mohr's salt solution = 500 ml

Therefore, amount of Fe $^{+2}$ in 10 ml of Mohr's salt solution = 0.05584 × V_{KMnO4} × S_{KMnO4} gm Amount of Fe $^{+2}$ in 500 ml of Mohr's salt solution = 0.05584 × V_{K2Cr2O7} × S_{K2Cr2O7} × 50 gm = 0.05584 × 7.5 × 0.06 × 50 gm = 1.2564 gm (Ans)

28. Calculate the known value of iron when 2g Mohr's salt is dissolved in 100 ml of solution. (Atomic weight of Fe = 55.84, Mol. weight of FeSO₄.(NH₄)₂SO₄.6H₂O = 392.14). Ans: 0.2847 gm

Solution: We know, $\frac{At.wt.of\ Cu\ /\ Fe \times Amount\ of\ salt\ taken\ (in\ gm)}{Mol.wt.of\ blue-vitriol\ /\ Mohr's\ salt}$

Given that, Atomic weight of Fe = 55.84, Mol. weight of FeSO₄. $(NH_4)_2SO_4.6H_2O = 392.14$, Amount of Mohr's salt taken = 2 gm;

Therefore, Known value of iron in 100 ml of blue vitriol solution = $(55.84 \times 2)/392.14$ = 0.28479 gm (Ans)