




EXPS: Determination of strength of a weak acid (CH_3COOH) against a strong alkali (NaOH) solution by measuring conductance

Name: Nafinur Leo, Id: 20-42195-1, Section: U

Theory:

Method involved: Conductometric titration

Reactions:



Here, my serial number 10

Experimental Data:

Table: Conductance-measurement of CH_3COOH and NaOH solutions using conductivity meter.

No. of reading	Vol. of CH_3COOH (in ml)	Vol. of NaOH (burette reading) (in ml)				Conductance (μS), γ
		Initial	Final	Difference	Total \times	
1	10	0	0	0	0	102
2	10	0	2	2	2	96
3	10	2	4	2	4	117
4	10	4	6	2	6	144
5	10	6	8	2	8	166
6	10	8	10	2	10	188
7	10	10	11	1	11	212
8	10	11	12	1	12	260
9	10	12	13	1	13	321
10	10	13	14	1	14	384
11	10	14	15	1	15	466
12	10	15	16	1	16	550

Calculations:

The strength of supplied CH_3COOH solution;

$$N_{\text{CH}_3\text{COOH}} \times V_{\text{CH}_3\text{COOH}} = N_{\text{NaOH}} \times V_{\text{NaOH}}$$

10 ml

0.1 N


$$\Rightarrow N_{\text{CH}_3\text{COOH}} = \frac{10 \times 0.1}{10}$$

$$= 0.1 \text{ N}$$

Results: The strength of supplied CH_3COOH

solution is 0.1 N

Conductometric Titration Curve CH_3COOH vs NaOH



Roll No.